

TUNNELWORKS

A2 CHEMISTRY

STUDENT SHEET 2

Oxygen detectives – Winkler equations

Take a sample of water from the river then:

Stage 1: Add manganese sulphate and then potassium iodide in potassium hydroxide.

First, the manganese sulphate reacts with the potassium hydroxide to form manganese hydroxide:



The manganese hydroxide reacts with the dissolved oxygen to form manganese oxyhydroxide.



(1 mol oxygen causes 2 mol manganic oxide to form)

The manganese was oxidised. Nothing's happened to the potassium iodide yet.

.....

Stage 2: Add sulphuric acid.

This dissolves the manganic oxide to form manganic sulphate, which reacts with the potassium iodide:



(2 mol manganic oxide releases 2 mol I₂ hence 2 mol I₂ is equivalent to 1 mol oxygen)

The manganese has been reduced again.

.....

Stage 3: Titrate with sodium thiosulphate using starch indicator.



(2 mol I₂ requires 4 mol thiosulphate. Hence 4 mol thiosulphate equates to 1 mol oxygen.)

.....

Stage 4: Quantify the results.

As an example, using 0.025 M thiosulphate solution, how much dissolved O₂ does each ml (cm³) of thiosulphate represent in the sample?

$1/1000 \times 0.025 \text{ M} \times 1/4 = 0.000625 \text{ mol O}_2 \text{ per ml } 0.025\text{M thiosulphate.}$

$\text{O}_2 = 32\text{g per mol so } 0.000625 \text{ mol} = 0.0002\text{g, or } 0.2\text{mg O}_2 \text{ per ml thiosulphate.}$

Scale the sample volume up to 1 litre to obtain a standard measurement of mg per litre.